



Recombinant *Saccharomyces cerevisiae* Poly (A) ribonuclease POP2

Product Code	CSB-EP330965SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39008
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MQSMNVQPRV LAVGGEQFFS QRQASEQHQQ QNMGPQVYSP KVNRRARMFPQ GMPVNTINGS VNQEMNNAYL LKQKNEPLLT QQQQQQQQQQ QPFNIGTPVS VASLPPGLNV LQQQQQQQQQ QQQQQQGVGL NRPLASQLPK HLTNQSMPPI FLPPPNYLFV RDVWKSNLYS EFAVIRQLVS QYNHVSISTE FVGTLARPIG TFRSKVDYHY QTMRANVDFL NPIQLGLSLS DANGNKPDNG PSTWQFNFEF DPKKEIMSTE SLELLRKSGI NFEKHENLGI DVFEFSQLLM DSGLMMDDSV TWITYHAAYD LGFLINILMN DSMPNNKEDF EWWWVHQYMPN FYDLNLVYKI IQEFKNPQLQ QSSQQQQQQQ YSLTTLADEL GLPRFSIFTT TGGQSLMLL SFCQLSKLSM HKFPNGTDFA KYQGVIYGID GDQ
Source	E.coli
Target Names	POP2
Protein Names	Recommended name: Poly(A) ribonuclease POP2 EC= 3.1.13.4 Alternative name(s): CCR4-associated factor 1
Expression Region	1-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.