



Recombinant Enterobacteria phage T4 Uncharacterized 22.0 kDa protein in e-segB intergenic region (y06Q)

Product Code	CSB-BP330991EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39224
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MNTLKKIVEF IRTKLGSAAMA KNLSVEEQYT AAAAKLLDKI KDLKTASVKS INEEKRIREL IVEKNKQAES KEREIRKLLS EGQDVTMHAK LGLLYRRTAE QLTTKADGYA EMRIEIAKKV VELDDARQEL AVKLEYIRET RAANALGIST ADDVVEIAAL TKVDIEDTLA RVETFNGNIS GVETTSADVQ EYINSLK
Source	Baculovirus
Target Names	y06Q
Protein Names	Recommended name: Uncharacterized 22.0 kDa protein in e-segB intergenic region
Expression Region	1-197
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.