



Recombinant Enterobacteria phage T4 Uncharacterized 39.0 kDa protein in nrdC-mobD intergenic region (y04P)

Product Code	CSB-MP330998EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39255
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MTRNEYIKSF NSVIDDKAIP MFGQNSVLSI INQWLNSVDA SIVSSTKFIH EIRKISSRVD KDVIKKTFKE SRLLSYLVNR DILGNFGKEI KRTKDVGYN WFGDVNSYHL NNKEDPENIF TRRWISNFR LFKKQILKSAS KLCYGDYRQI HPLASDMIII KEYELDKNKV SIFVNYGFFT PETNQKNINK FFSIASTITR QLETALLCME TVENIHTYPF KNICGWEGYK LVISLREVKC AYSPTSKEIY QQKCDEIVNT PKEETTLEEL MECLDDSP EP IEIRPEVIAL EKAYKEVLEI SNKAQKEYEQ AKKIWEESVN RLDRLEQALQ LIK
Source	Mammalian cell
Target Names	y04P
Protein Names	Recommended name: Uncharacterized 39.0 kDa protein in nrdC-mobD intergenic region
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.