



Recombinant Enterobacteria phage T4 NAD--protein ADP-ribosyltransferase modA (modA)

Product Code	CSB-BP331019EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39421
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MKYSVMQLKD FKIKSMDASV RASIREELLS EGFNLSEIEL LIHCITNKPD DHSWLNEIHK SRLVPNDKPL WRGVPAETKQ VLNQGIDIIT FDKVVSASYD KNIALHFASG LEYNTQVIFE FKAPMVFNFQ EYAIKALRCK EYNPNFKFPD SHRYRNMELV SDEQEVMIPIA GSVFRIADRY EYKKCSTYTI YTLDFEGFNL
Source	Baculovirus
Target Names	modA
Protein Names	Recommended name: NAD--protein ADP-ribosyltransferase modB EC= 2.4.2.-
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.