



Recombinant *Saccharomyces cerevisiae* 3-oxoacyl-[acyl-carrier-protein] synthase homolog (CEM1)

Product Code	CSB-EP331037SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39525
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSRRVVITGL GCVTPLGRSL SESWGNLLSS KNGLTPITSL PNYNEDYKLR EKSIPSTITV GKIPENFQNE NSAINKLLFT SQDERRTSSF IKLALRTTYE ALHNAGLLNP NDITINTSLC NLDHFGCLIG SGIGSIQDIY QTSLQFHNDN KRINPYFVPK ILTNMAAGNV SIKFNLRGLS HSVSTACATG NNSIGDAFNF IRLGMQDICV AGASETSLHP LSLAGFIRAK SITTINGISRP FDTQRSGFVL GEGCGMIVME SLEHAQKRNA NIISELVGYG LSSDACHITS PPADGNGAKR AIEMALKMAR LEPTDVDYVN AHATSTLLGD KAECLAVASA LLPGRSKSKP LYISSNKGAI GHLLGAAGAV ESIFTICSLK DDKMPHTLNL DNVLTLENNE ADKLFHIRDK PIVGANPKYA LCNSFGFGGV NTSLLFKKWE GS
Source	E.coli
Target Names	CEM1
Protein Names	Recommended name: 3-oxoacyl-[acyl-carrier-protein] synthase homolog EC=2.3.1.41 Alternative name(s): Beta-ketoacyl-ACP synthase homolog
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.