



Recombinant Escherichia coli D-galactarate dehydratase (garD)

Product Code	CSB-BP331084ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39829
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MANIEIRQET PTAFYIKVHD TDNVAIIVND NGLKAGTRFP DGLELIEHIP QGHKVALLDI PANGEIIRYG EVIGYAVRAI PRGSWIDESM VVLPEAPPLH TLPLATKVPE PLPPLEGYTF EGYRNADGSV GTKNLLGITT SVHCVAGVVD YVVKIIERDL LPKYPNVDGV VGLNHLYGCG VAINAPAAVV PIRTIHNISL NPNFGGEMV IGLGCEKLQP ERLLTGTDDV QAIPVESASI VSLQDEKHVG FQSMVEDILQ IAERHLQKLN QRQRETCPAS ELVVMQCGG SDAFSGVTAN PAVGYASDLL VRCGATVMFS EVTEVRDAIH LLTPRAVNEE VGKRLEEME WYDNYLNMGK TDRSANPSPG NKKGGLANVV EKALGSIKAKS GKSAIVEVLS PGQRPTKRGL IYAATPASDF VCGTQQVASG ITVQVFTTGR GTPYGLMAVP VIKMATRTEL ANRWFDLMDI NAGTIATGEE TIEEVGWKLF HFILDVASGK KKTFSQWGL HNQLAVFNPA PVT
Source	Baculovirus
Target Names	garD
Protein Names	Recommended name: D-galactarate dehydratase Short name= GalcD EC= 4.2.1.42
Expression Region	1-523
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.