



# Recombinant *Saccharomyces cerevisiae* Putative uncharacterized protein YER076C (YER076C)

<b>Product Code</b>	CSB-EP331119SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P40049
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	EIYLTGG ESQGLPSGFV QMDDDLAIAP VSMVEFYQTI GLTANGTVPE SFNKRDATEY PNIISNITTQ AANFTQHSLV EQLQNDVTAI SISNAINVAV DGSVETPADL QYNRNQETGE STLCRAKFGY VEVRTWSRLY NRAVSSTTLT TNLNSYIAQW VAWAVHGSGD KKFCGSQEFT NIFFDGQEGW SLFVKTWSTN SSCDITASEG NLTC AVRVS V SSMH NHGKTA FCV TYSHGDS WRAELRVVAN DAWSHYYPWS IDCPEVDKNN MAINDCFDQA QG
<b>Source</b>	E.coli
<b>Target Names</b>	YER076C
<b>Protein Names</b>	Recommended name: Putative uncharacterized protein YER076C
<b>Expression Region</b>	24-302
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.