



Recombinant *Saccharomyces cerevisiae* Putative uncharacterized protein YER076C (YER076C)

Product Code	CSB-MP331119SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40049
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	EIYLTGG ESQGLPSGFW QMDDDLAIAP VSMVEFYQTI GLTANGTVPE SFNKRDATEY PNIISNITTQ AANFTQHSLV EQLQNDVTAI SISNAINVAV DGSVETPADL QYNRNQETGE STLCRAKFGY VEVRTWSRLY NRAVSSTTLT TNLNSYIAQW VAWAVHGSGD KKFCGSQEFT NIFFDGQEGW SLFVKTWSTN SSCDITASEG NLTC AVRVS SV SSMHNHGKTA FCVTYSHGDS WRAELRVVAN DAWSHYYPWS IDCPEVDKNN MAINDCFDQA QG
Source	Mammalian cell
Target Names	YER076C
Protein Names	Recommended name: Putative uncharacterized protein YER076C
Expression Region	24-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.