



Recombinant *Saccharomyces cerevisiae* Protein MMF1, mitochondrial (MMF1)

Product Code	CSB-YP331139SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40185
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	ITT LTPVSTKLAP PAAASYSQAM KANNFVYVSG QIPYTPDNKP VQGSISEKAE QVFQNVKNIL AESNSSLDNI VKVNVFLADM KNFAEFNSVY AKHFHTHKPA RSCVGVASLP LNVDLEMEVI AVEKN
Source	Yeast
Target Names	MMF1
Protein Names	Recommended name: Protein MMF1, mitochondrial Alternative name(s): Isoleucine biosynthesis and maintenance of intact mitochondria 1 Maintenance of mitochondrial function 1
Expression Region	18-145
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.