



Recombinant Human Wilms tumor protein 1-interacting protein (WTIP)

Product Code	CSB-EP026161HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NIX2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQRSRAGADE AALLLAGLAL RELEPGCGSP GRGRRGPRPG PGDEAAPALG RRGKGS GGPE AGADGLSRGE RGPRRAAVPE LSAQPAGSPR ASLAGSDGGG GGSARSSGI SLGYDQRHGS PRSGRSDPRP GPGPPSVGSA RSSVSSLGSR GSAGAYADFL PPGACPAPAR SPEPAGPAPF PLPALPLPPG REGGPSAAER RLEALTRELE RALEARTARD YFGICIKCGL GIYGAQQACQ AMGSLYHTDC FTCDSCGRRRL RGKAFYNVGE KVYCQEDFLY SGFQQTADKC SVCGLHIMEM ILQALGKSYH PGCFRCSVCN ECLDGVFPFTV DVENNIYCVR DYHTVFAPKC ASCARPILPA QGCETTIRVV SMDRDYHVAC YHCEDCGLQL SGEEGRRCYP LAGHLLCRRC HLRRLQPGPL PSPTVHVTEL</p>
Source	E.coli
Target Names	WTIP
Protein Names	Recommended name: Wilms tumor protein 1-interacting protein Short name= WT1-interacting protein
Expression Region	1-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.