



# Recombinant Human Wilms tumor protein 1-interacting protein (WTIP)

<b>Product Code</b>	CSB-MP026161HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A6NIX2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQRSRAGADE AALLLAGLAL RELEPGCGSP GRGRRGPRPG PGDEAAPALG RRGKGS GGPE AGADGLSRGE RGPRAAVPE LSAQPAGSPR ASLAGSDGGG GGSARSSGI SLGYDQRHGS PRSGRSDPRP GPGPPSVGSA RSSVSSLGSR GSAGAYADFL PPGACPAPAR SPEPAGPAPF PLPALPLPPG REGGPSAAER RLEALTRELE RALEARTARD YFGICIKCGL GIYGAQQACQ AMGSLYHTDC FTCDSCGRRRL RGKAFYNVGE KVYCQEDFLY SGFQQTADKC SVCGLHIMEM ILQALGKSYH PGCFRCSVCN ECLDGVPFV DVENNIYCVR DYHTVFAPKC ASCARPILPA QGCETTIRVV SMDRDYHVAC YHCEDCGLQL SGEEGRRCYP LAGHLLCRRC HLRRLQGPL PSPTVHVTEL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	WTIP
<b>Protein Names</b>	Recommended name: Wilms tumor protein 1-interacting protein Short name= WT1-interacting protein
<b>Expression Region</b>	1-430
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.