



Recombinant Human Proline-rich protein 19 (PRR19)

Product Code	CSB-BP018792HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NJB7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDTQGPVSQLP FQQPEKPGRV RRRKTRRERN KALVGSRRPL AHHDPPVAIR DPPVVPTASK LVVITQGRLS REHRGLFNHE VKSLDVARLL SSGTLVPGSP TLPKPSPPS GRAQEPAPRS RDKENQVPGG SGPGPSPE LSGVGQLLAE LQCQLSLPQA FRRNLIQDA RDAIVHTLQA CHGCVPDLAL VLRGCQPPLP GAKPGVSEK MTPFWINSPD QVPEQERQRK QQGTFEFTFP MPYTSSMPTA HRGSLAPPRG PWPPYFPSLS SPSGTAWGPP TAFDLLKSIW LVATPPPPRP WGVGLPQPLP QPSSPLLPRV SVLDWSPSPS SPLPSLSWVV AQSSPEAWSF PPMRLY
Source	Baculovirus
Target Names	PRR19
Protein Names	Recommended name: Proline-rich protein 19
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.