



Recombinant *Saccharomyces cerevisiae* Cop9 signalosome complex subunit 11 (PCI8)

Product Code	CSB-EP331184SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40512
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MFHGAKGPLL IERIGHLLSI NYGEKEERKR WAMQGISYLQ EVQCTSTPYL EILVEESGLR PVSLLNSQLV GKPHFSLG FDENIARDII SHNFQNAIFQ MESEEVPLTK RYQHLEKITQ ISLLCKNFKG IEEIEYNVKN IIQGRKNFDM LNSMEKDRIS HEVVQDDSFSL LRIQMMLLCV SYFLQERYFD CTKFFTTMMT SEPLTLKVLS EHLDCMNFIS KEEFIMMVNI SVLISIPLDN YDDFIYLSDL KQFFQMTPLL VNCLELLINT NFNKFFKIWH GEINKICMES LFLEPSWSSS AAVIMRCKIY FFYLRISKKL QFSYLSSTLG IDLEDIKEEL TKLIISGQLN FEIDGDVIHF EDSSILQSIV NEISRNGTMI NEVIDKLNKNE NTDLKDIIQG NPLMYSGGNN TATIINNESS DDMDIDEVND RSDISDSEGG LFEC
Source	E.coli
Target Names	PCI8
Protein Names	Recommended name: Cop9 signalosome complex subunit 11 Alternative name(s): Proteasome-COP9 signalosome-eIF3 domain protein 8
Expression Region	1-444
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.