



Recombinant *Saccharomyces cerevisiae* Zinc finger protein RTS2 (RTS2)

Product Code	CSB-MP331248SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40962
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MADYDSAKYW SKQGARRGLQ KTRYQCICQ RQCKDANGFQ SHNKSPSHLR KISQVTAEDA RRYNIQFEKG FLQLLKQRHG EKWIDANKVY NEYVQDRDHV HMNATMHRSL TQFVRYLGRA GKVDVDMID DTSENVEGPL LIRIHPSSLS SPSEDGMLRS QQEEQEVIAA ELLKRQLNRA KRQTEKVYQP EMKSEISGDS TLKRVQVTFH GNDRVNKKKK KVPPRKDGKIK FR
Source	Mammalian cell
Target Names	RTS2
Protein Names	Recommended name: Zinc finger protein RTS2
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.