



Recombinant Arabidopsis thaliana 60S ribosomal protein L13-1 (RPL13B)

Product Code	CSB-BP331266DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P41127
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MKHNNVIPNG HFKKHWQNYV KTWFNQPARK TRRRIARQKK AVKIFPRPTS GPLRPVVHGQ TLKYNMKVRT GKGFTLEELK AAGIPKKLAP TIGIAVDHRR KNRSLEGLQT NVQRLKTYKT KLVIFPRRAR KVKAGDSTPE ELANATQVQG DYLPIVREKP TMELVKLTSE MKSFKAFDKI RLERTNKRHA GARAKRAAEA EKEEEK
Source	Baculovirus
Target Names	RPL13B
Protein Names	Recommended name: 60S ribosomal protein L13-1 Alternative name(s): Protein BBC1 homolog
Expression Region	1-206
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.