



Recombinant *Saccharomyces cerevisiae* Target of rapamycin complex subunit LST8 (LST8)

Product Code	CSB-EP331278SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P41318
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSVILVSAGY DHTIRFWEAL TGVCSTRTIQH SDSQVNRLEI TNDKLLATA GHQNVRLYDI RTTNPNPVAS FEGHRGNVTS VSFQQDNRWM VTSSDGTIK VWDVRSPSIP RNYKHNPVN EVVIHPNQGE LISCDRDGNI RIWDLGENQC THQLTPEDDT SLQSLSMASD GSMLAAANTK GNCYVWEMPEN HTDASHLKPV TKFRAHSTYI TRILLSSDVK HLCATCSADHT ARVWSIDDDF KLETTLDGHQ RWVWDCAFSA DSAYLVTASS DHYVRLWDL TREIVRQYGG HHKGAVCVL NDV
Source	E.coli
Target Names	LST8
Protein Names	Recommended name: Target of rapamycin complex subunit LST8 Short name= TORC subunit LST8 Alternative name(s): Lethal with SEC13 protein 8
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.