



Recombinant *Saccharomyces cerevisiae* NADPH dehydrogenase 3 (OYE3)

Product Code	CSB-MP331351SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P41816
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPFVKGFEP I SLRDTNLFEP IKIGNTQLAH RAVMPPLTRM RATHPGNIPN KEWAAVYYGQ RAQRPGTMII TEGTFISPQA GGYDNAPGIW SDEQVAEWKN IFLAIHDCQS FAWVQLWSLG WASFPDVLAR DGLRYDCASD RVYMNATLQE KAKDANNLEH SLTKDDIKQY IKDYIHAANK SIAAGADGVE IHSANGYLLN QFLDPHSNKR TDEYGGTIEN RARFTLEVVD ALIETIGPER VGLRLSPYGT FNSMSGGAEP GIIAQYSYVL GELEKRAKAG KRLAFVHLVE PRVTDPSLVE GEGEYSEG TN DFAYSIWKGP IIRAGNYALH PEVVREQVKD PRTLIGYGRF FISNPDLVYR LEEGLPLNKY DRSTFYTMSA EGYTDYPTYE EAVDLGWKNK
Source	Mammalian cell
Target Names	OYE3
Protein Names	Recommended name: NADPH dehydrogenase 3 EC= 1.6.99.1 Alternative name(s): Old yellow enzyme 3
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.