



# Recombinant Schizosaccharomyces pombe 26S protease regulatory subunit 8 homolog (let1)

<b>Product Code</b>	CSB-EP331354SXV-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P41836
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTEVLKTNVL QSNENIVQYY TQKIQDAELA ILQKTQNLRR LEAQRNGLNA RVRLLRREEIQ LLQEPGSYVG EVIKTMGKNK VLVKVHPEGK YVVDISPDID IKEIKPNIRV ALRNDYQLI KILPNKVDPL VSLMMVEKIP DSTYEMVGGL EKQIKEIKEV IELPVKHPEL FESLGIPQPK GILLYGPPGT GKTLLARAVA HHTDCKFIRV SGSELVQKYI GEGSRMVREL FVMAREHAPS IIFMDEIDSI GSSRSOSSGG SGDSEVQRTM LELLNQLDGF EATKNIKVIM ATNRIDILDP ALLRPGRIDR KIEFPPPSAE ARAEILRIHS RSMNLTRGID LKSIAEKMNG ASGAELKGVC TEAGMFALRE RRVHVTQEDF ELAVAKVLNK GDSGEMSLQK LFK
<b>Source</b>	E.coli
<b>Target Names</b>	let1
<b>Protein Names</b>	Recommended name: 26S protease regulatory subunit 8 homolog Alternative name(s): Protein let1
<b>Expression Region</b>	1-403
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.