



# Recombinant *Saccharomyces cerevisiae* Mannose-1-phosphate guanyltransferase (PSA1)

<b>Product Code</b>	CSB-BP331374SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P41940
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKGLILVGGY GTRLRPLTLT VPKPLVEFGN RPMILHQIEA LANAGVTDIV LAVNYRPEVM VETLKKYEKE YGVNITFSVE TEPLGTAGPL KLAEDVLKGD NSPFFVLNSD VICEYPFKEL ADFHKAHGGK GTIVATKVDE PSKYGVIVHD IATPNLIDRF VEKPKFVGN RINAGLYILN PEVIDLIEMK PTSIEKETFP ILVEEKQLYS FDLEGFWMDV GQPKDFLSGT VLYLNSLAKR QPKKLATGAN IVGNALIDPT AKISSTAKIG PDVVIGPNVT IGDGVRITRS VVLCNSTIKN HSLVKSTIVG WNSTVGQWCR LEGVTVLGDD VEVKDEIYIN GSKVLPKHSI SDNVPKEAII M
<b>Source</b>	Baculovirus
<b>Target Names</b>	PSA1
<b>Protein Names</b>	Recommended name: Mannose-1-phosphate guanyltransferase EC= 2.7.7.13 Alternative name(s): ATP-mannose-1-phosphate guanylyltransferase GDP-mannose pyrophosphorylase NDP-hexose pyrophosphorylase
<b>Expression Region</b>	1-361
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.