



Recombinant Caenorhabditis briggsae Apoptosis regulator ced-9 (ced-9)

Product Code	CSB-MP331378CXX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P41957
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	MVDSMDMANS SQNTFRRRTM ATSEMREFLS TKDAEPNNFG MQTIPESPTP STPTRRMSIG DSTRIYDWEE PRFNIQGFVV DYFTYRIAQN GLDWYDAPAL PDGVQKEHEM MRSLGTIFEK RHMEMFENFS EQLLAVPKIS FSLYQEVVQT VGNSSNTPCP MSYGRILIGLI SFGGMVAAKM MESAELQGQV RNLLMYTSLF IKTRIRQSWK EHNRSWADFM KLGQQMKEDY EKEKDAEEGK RLKSWSIIGA SVIAVIVCGR IIFSFK
Source	Mammalian cell
Target Names	ced-9
Protein Names	Recommended name: Apoptosis regulator ced-9 Alternative name(s): Cell death protein 9
Expression Region	1-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.