



Recombinant Putative casein kinase I C03C10.2 (C03C10.2)

Product Code	CSB-MP331433CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42169
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MNKSKEAETD ITQIFPNLDI LPKKFSAKEL AERRIKKLHV APGSIFMNRW SIEGVINGGG YGQIFMVM DV KKNDERAMKI EPKLRAEVIT KRMIMEQQVL MKMQGKTHIP TMYASGFNDQ FNFIIMQLLS MNVGDFRKRS PLGRLSKETV GRIAYQTLNA LKDIHDMGYV HRDVK PANIC FGVHAQNRHI LYLLDFGLVR RFKTESGVC I PWRINAGFKG TERYVSVRVH EKLEQTPWDD AFSVLYTAYE LVVGELPWRY LEDIHEIHGV KKL MNEVTKN GEMFKDIASI LVDFHKMILE CDPVVELPYE K LLECLKCLY SPKS LLEPYD WEDGYKTTLN I
Source	Mammalian cell
Target Names	C03C10.2
Protein Names	Recommended name: Putative casein kinase I C03C10.2 EC= 2.7.11.1
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.