



Recombinant Uncharacterized protein C03C10.4 (C03C10.4)

Product Code	CSB-EP331434CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42171
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MTEGDDQLIS IRKQLENLNN ATDDINSYEM KLETVKKQFC ETQLMFNKEM LGIPKKLAKH ISKSRQFFDL KSRESEIRRC VQQAAAQFER QKTSVEMARE QVQILHNSLN NNQELDAEKQ YVDVIEQQLE LVKEAEGECL KAEKCHASRV RDLLQLEMAL RKCLEENGSA IKKSRPYER KEVLTRTMNS QLELMSILEH EVQERKDSYS DSMRALEQIS DQIHQERSSQ SSLAPSSDAE SDSS
Source	E.coli
Target Names	rei-1
Protein Names	Recommended name: Uncharacterized protein C03C10.4
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.