



# Recombinant Capsule polysaccharide export protein kpsC (kpsC)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP331454ENL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Uniprot No.</b>       | P42217                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Immunogen Species</b> | Escherichia coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Sequence</b>          | MIGIYSPGIW RIPHLEKFLA QPCQKLSLLR PVPQEVDAIA VWGHRPSAAK<br>PVAIAKAAGK PVIRLEDGFV RSLDLGVNGE PPLSLVDDY CIYYDASKPS<br>ALEKLVQDKA GNTALISQAR EAMHTIVTGD LSKYNLAPAF VADESERSDI<br>VLVVDQTFND MSVTYGNAGP HEFAAMLEAA MAENPQAEIW VKVHPDVLEG<br>KKTGYFAALR ATQVRVLIIE NVSPQSLLRH VSRVYVVTSS YGFEALLAGK<br>PVTCFGQPWY AGWGLTDDRH PQSALLSARR GSATLEELFA AAYLRYCRYI<br>DPQTGEVSNL FTVLQWLQLQ RRHLQQRNGY LWAPGLTLWK SAILKPFLQT<br>ATNRLSFSRR CTAASACVWV GVKGEQQWRA EAQRKSLPLW<br>RMEDGFLRSS GLGSDLLPPL SLVLDKRGYI YDATRPSDLE VLLNHSQTL<br>AQQMRAEKLK QRLVESKLSK YNLGADFSLP AEAKDKKIIL VPGQVEDDAS<br>IKTGTVSIKS NLELLRTVRE RNPWAYIYK PHPDVLVGNR KGNIPTELIA<br>ELADYQALDA DIIQCIQRAD EVHTMTSLSG FEALLHGKQV HCYGLPFYAG<br>WGLTVDEHHC PRREQKLTIA DLIYQTLIVY PTYIHPTLLQ PISVEEAAEY<br>LIQTPRKPMF ITRKKAGRVI RYYRKLIMFC KVRFG |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Target Names</b>      | kpsC                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Names</b>     | Recommended name: Capsule polysaccharide export protein kpsC                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Expression Region</b> | 1-675                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.