



Recombinant Bacillus subtilis Protein ioIH (ioIH)

Product Code	CSB-YP331529BRJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42418
Product Type	Recombinant Protein
Immunogen Species	Bacillus subtilis (strain 168)
Purity	≥85% (SDS-PAGE)
Sequence	MKLALDPSMY RDDLTLEEMV YKTAELGYEY IELSPREDFC PFYKYPKVDS AKIKQLKRL RDTGVKLSSL LPLYHWAGPD EDRRQAAVRN WKRAIEIAVE LEV DLMNSEF SGSKYDPLTS EEKFIKSMDE LLPVFEKEGI KLN LQAHPYD FIETHKGAMD MIRALDKDWI NLVYSTAH TF FYDDGKGDIE TMFDEAGDRL THVLFADTYN HKA AHGLRYI VNPPDAKVTV HQHLDIGQQGE VDFDTIFRKL REMKFDGIAT NAVFAWVDER ADESSRFMLK KMKEELGLA
Source	Yeast
Target Names	ioIH
Protein Names	Recommended name: Protein ioIH
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.