



# Recombinant Human Zinc finger protein 98 (ZNF98)

<b>Product Code</b>	CSB-MP027106HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A6NK75
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPGPLGSLEM GVLTRFDVAL EFSLEEWQCL DTAQQNLYRN VMLENYRNLV FVGIAASKPD LITCLEQGKE PWNVKRHEMV TEPPVVYSYF AQDLWPKQGK KNYFQKVILR TYKKCGRENL QLRKYCKSMD ECKVHKECYN GLNQCLTTTQ NKIFQYDKYV KVFHKFSNSN RHKIGHTGKK SFKCKECEKS FCMLSHLAQH KRIHSGEKPY KCKEKGKAYN EASNLSTHKR IHTGKKPYKC EECGKAFNRL SHLTTHKIIH TGKKPYKCEE CGKAFNQSAN LTTHKRIHTG EKPYKCEECEG RAFSQSSTLT AHKIIHAGEK PYKCEECEGKA FSQSSTLTTH KIIHTGEKPY KCEECEGKA FS RLSHLTTHKR IHSGEKPYKC EECGKAFKQS STLTTHKRIH AGEKPYKCEV CSKAFSRFSH LTTHKRIHTG EKPYKCEECEG KAFNLSSQLT THKIIHTGKPYKCEECEGKA FNQSSTLSKH KVIHTGKPYKCEECEGKAFN QSSHLTTHKM IHTGKPYKC EECGKAFNNS SILNRHKMIH TGEKLYKPES CNNACDNI AK ISKYKRNCAG EK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF98
<b>Protein Names</b>	Recommended name: Zinc finger protein 98 Alternative name(s): Zinc finger protein 739 Zinc finger protein F7175
<b>Expression Region</b>	1-572
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.