



Recombinant Arabidopsis thaliana G-box-binding factor 3 (GBF3)

Product Code	CSB-EP331632DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42776
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGNSSEPKP PTKSDKPSSP PVDQTNVHVY PDWAAMQAYY GPRVAMPPYY NSAMAASGHP PPPYMWNPQH MMSPYGAPYA AVYPHGGGVY AHPGIPMGSL PQQQKDPPLT TPGTLLSIDT PTKSTGNTDN GLMKKLKEFD GLAMSLGNGN PENGADEHKR SRNSSETDGS TDGSDGNTTG ADEPKLKR SR EGTPTKDGKQ LVQASSFHSV SPSSGDTGVK LIQSGAILS PGVSANSNPF MSQSLAMVPP ETWLQNEREL KRERRKQSNR ESARRSRLRK QAETEELARK VEALTAENMA LRSELNQLNE KSDKLRGANA TLLDKLKCSE PEKRV PANML SRVKNSGAGD KNKNQGDND S NSTSKLHQLL DTKPRAKAVA AG
Source	E.coli
Target Names	GBF3
Protein Names	Recommended name: G-box-binding factor 3 Alternative name(s): bZIP transcription factor 55 Short name= AtbZIP55
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.