



Recombinant Escherichia coli Uncharacterized fimbrial-like protein yraH (yraH)

Product Code	CSB-MP331681ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42913
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	TDGEIVFD GEILKSACEI NDSDKKIEVA LGHYNAEQFR NIGERSPKIP FTIPLVNCMP TGWEHDNGNV EASFRLWLET RDNGTVPNFP NLAKVGSFAG IAATGVGIRI DDAESGNIMP LNAMGNDNTV YQIPAESNGI VNVDLIAYYV STVVPSEITP GEADAI NVVT LDYR
Source	Mammalian cell
Target Names	yraH
Protein Names	Recommended name: Uncharacterized fimbrial-like protein yraH
Expression Region	23-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.