



# Recombinant *Saccharomyces cerevisiae* Uncharacterized GTP-binding protein YGR210C (YGR210C)

<b>Product Code</b>	CSB-MP331692SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P42942
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPRDPLIGIV GKPSSGKSTT LNSLTDAGAA VGAFPFTTIE PNQATGYLQV ECACSRFGKE DLCKPNYGWC SKGKRHIPIK LLDVAGLVPG AHSGRGLGNK FLDDL RHADA LIHVVDVSGT TDAEGKNTRG YDPLNDIEWL QDEIRLWVEG NLKKRWGSIV RKHTATKSSI VDTLQAQFGG YGSHAPMIQK ALDRLKGLPP LEKWDDWIT RVVKSFMMEK FPTVLALNKI DHPDADKNVS KIMLKYPDTK AVLTSAVTEV FLRKLKKQGF ILYEEGTEFV DTFEDEPKL KPLDDKILNR IENIRDLVLY RFGSTGVVQV LQAATDILGL IPVYTVKNIQ TFTGGNGTNV FRDCFLVKRG TPVGKVARYI MGGEVTIASI ETVGGVRVSE ESLVEPGKND ILGFKIAPRS A
<b>Source</b>	Mammalian cell
<b>Target Names</b>	YGR210C
<b>Protein Names</b>	Recommended name: Uncharacterized GTP-binding protein YGR210C
<b>Expression Region</b>	1-411
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.