



Recombinant Phleum pratense Pollen allergen Phl p 1 (PHLPI)

Product Code	CSB-BP331772EUQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43213
Product Type	Recombinant Protein
Immunogen Species	Phleum pratense (Common timothy)
Purity	>85% (SDS-PAGE)
Sequence	IPKVPPG PNITATYGDK WLDAKSTWYG KPTGAGPKDN GGACGYKDVD KPPFSGMTGC GNTPIFKSGR GCGSCFEIKC TKPEACSGEP VVVHITDDNE EPIAPYHFDL SGHAFGAMAK KGDEQKLRSA GELELQFRRV KCKYPEGTKV TFHVEKGSNP NYLALLVKYV NGDGDVVAVD IKEKGKDKWI ELKESWGAIW RIDTPDKLTG PFTVRYTTEG GTKTEAEDVI PEGWKADTSY ESK
Source	Baculovirus
Target Names	PHLPI
Protein Names	Recommended name: Pollen allergen Phl p 1 Alternative name(s): Allergen Phl p I Allergen= Phl p 1
Expression Region	24-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.