



Recombinant Arachis hypogaea Allergen Ara h 1, clone P17, partial

Product Code	CSB-MP331779ANE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43237
Product Type	Recombinant Protein
Immunogen Species	Arachis hypogaea (Peanut)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	KSPYRK TENPCAQRCLQSCQQEPDDLKQKACESRCKLEYDPRCVYDTGAT NQRHPPGERTGRQPGDYDDRRQPRREEGGRWGPAEPREREREEDWRQ PREDWRRPSHQQPRKIRPEGREGEQEWGTPGSEVREETSRRNNPFYFSSRRF STRYGNQNGRIRVLQRFQDQSKQFQNLQNHRIVQIEARPN
Source	Mammalian cell
Protein Names	Allergen Ara h I Allergen: Ara h 1
Expression Region	26-216aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.