



Recombinant Arabidopsis thaliana Probable cysteine proteinase A494 (At2g21430)

Product Code	CSB-EP331790DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43295
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	LPEEFDWRD RGAVTPVKNQ GSCGSCWSFS TTGALEGAHF LATGKLVSLSEQQLVDCDHE CDPEEEGSCD SGCNGLMNS AFEYTLKTGGLMREKDYPYT GTDGGCKLD RSKIVASVSN FSVVSINEDQ IAANLIKNGPLAVAINAAYM QTYIGGVSCP YICSRRLNHG VLLVGYGSAG FSQARLKEKPYWIIKNSWGE SWGENGFYKI CKGRNICGVD SLVSTVAATT S
Source	E.coli
Target Names	RD19B
Protein Names	Recommended name: Probable cysteine proteinase A494 EC= 3.4.22.-
Expression Region	132-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.