



Recombinant *Saccharomyces cerevisiae* Pre-mRNA-splicing factor SAD1 (SAD1)

Product Code	CSB-BP331880SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43589
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MEVDNKRKRRHS EDELKQEAVK KIKSQEPNYA YLETVVREKL DFDSEKICCI TLSPNLVYCC LVCGHYYQGR HEKSPAFIHS IDENHHVFLN LTSLKFYMLP QNVQILHDGE VQLLNSIKFA AYPTYCPKDL EDFPRQCDFL SNRTYLNFGI GFTNAATYDY AHSVLLISH MVPVRDHFLN NHFDNQGEFI KRLSICVKKI WSPKLFKHHL SVDDFVSYLK VREGLNLPNI DPRLFLLWLF NKICSSSNDL KSILNHSCKG KVKIAKVENK PEASESVTGG VIVKPFWVLT LDLPEFSPFE DGNSVDDL PQ INITKLLTKF TKSRSSSTST VFELTRLPQF LIFHFNRFDR NSDHPVKNRN QTLVEFSSEL EILHVKYRLK ANVVHVVIKQ PSTDGNAFNG DEKSHWITQL YDNKSEKWIE IDGINTTERE AELLFLKETF IQVWEKQE
Source	Baculovirus
Target Names	SAD1
Protein Names	Recommended name: Pre-mRNA-splicing factor SAD1 Alternative name(s): snRNP assembly-defective protein 1
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.