



Recombinant *Saccharomyces cerevisiae* Kinetochores-associated protein CNN1 (CNN1)

Product Code	CSB-BP331894SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43618
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSTPRKAAGN NENTEVSEIR TPFRRERALEE QRLKDEVLIR NTPGYRKLLS ASTKSHDILN KDPNEVRSFL QDLSQVLARK SQGNDTTTNK TQARNLIDEL AYEESQPEEN ELLRSRSEKL TDNNIGNETQ PDYTSLSQTV FAKLQERDKG LKSRKIDPII IQDVPTTGHE DELTVHSPDK ANSISMEVLR TSPSIGMDQV DEPPVRDPVP ISITQQEEPL SEDLPSSDKE ETEEAENEDY SFENTS DENL DDIGNDPIRL NVPARRSSI KPLQIMDLKH LTRQFLNENR IILPKQTWST IQEESLNIMD FLKQKIGTLQ KQELVDSFID MGIINNVDDM FELAHHELLPL ELQSRIESYL F
Source	Baculovirus
Target Names	CNN1
Protein Names	Recommended name: Kinetochores-associated protein CNN1 Alternative name(s): Co-purified with NNF1 protein 1
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.