



# Recombinant Human C2 calcium-dependent domain-containing protein 4B (C2CD4B)

<b>Product Code</b>	CSB-EP008080HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A6NLJ0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MRLLEKLCSS AAGSSAPKPA FAKVLTPNRI PEFCIPRLP APCTLESP AAVPRRCAA ESDLWPRAAD EDAGRTDWDW RSQAALSLPH LPRVRTTYGF CALLESPHTR RKESLLLGGP PAPERPAHSC GGGGGPDAPL GTLCGPRGPG PATPAAPGGP RLPQDALAAG PRRCRLLRVP DGLLSRALRA GRSRRLARVR SVSSGNEDEE RRAGSESPAR APSSSPLSSR APLPERLEAK GTVALGRAGD ALRLAAEYCP GTRRLRLRL RAESLFGGAP GPRAVRCRLS LVLRPPGTAR WQCSAVVGRS RKASFDQDFC FDGLSEDEV RLA VRVKARD EGRGRDRGRL LGQGELSLGA LLLL
<b>Source</b>	E.coli
<b>Target Names</b>	C2CD4B
<b>Protein Names</b>	Recommended name: C2 calcium-dependent domain-containing protein 4B Alternative name(s): Nuclear-localized factor 2 Protein FAM148B
<b>Expression Region</b>	1-364
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.