



Recombinant Mouse Dimethylaniline monooxygenase [N-oxide-forming] 3 (Fmo3)

Product Code	CSB-MP008749MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P97501
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	KKKVAIIGA GVSGLAAIRS CLEEGLEPTC FERSDDVGGL WKFSDHIEEG RASIYQSVFT NSSKEMMCFP DFPYPDDFPN FMHHSKLQEQY ITSFAKEKNL LKYIQFETPV TSINKCPNFS TTGKWEVTTE KHGKKETAVF DATMICSGHH IFPHVPKDSF PGLNRFKGGK FHSRDYKEPG IWKGKRVLVI GLGNSGCDIA AELSHVAQKV TISSRSGSWV MSRVWDDGYP WDMVVLTRFQ TFLKNNLPTA ISDWWYTRQM NARFKHENYG LVPLNRTLRLK EPVFNDLPA RILCGMVTIK PNVKEFTETS AVFEDGTMFE AIDCVIFATG YGYAYPFLDD SIKSRNNEV TLYKGVFPPQ LEKPTMAVIG LVQSLGATIP ITDLQARWAA QVIKGTCTLP SVNDMMDDID EKMGEKFKWY GNSTTIQTDY IVYMDELASF IGAKPNLLWL FLKDPRLAVE VFFGPCSPYQ FRLVGPGKWS GARNAILTQW DRSLKPMKTR VVSKVQKSCS HFYSRLLRLL AVPVLLIALF LVLI
Source	Mammalian cell
Target Names	Fmo3
Protein Names	Recommended name: Dimethylaniline monooxygenase [N-oxide-forming] 3 EC= 1.14.13.8 Alternative name(s): Dimethylaniline oxidase 3 Hepatic flavin-containing monooxygenase 3 Short name= FMO 3
Expression Region	2-534
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.