



Recombinant Cell division protein FtsZ (ftsZ)

Product Code	CSB-BP332608WAR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P45485
Product Type	Recombinant Protein
Immunogen Species	Wolbachia sp.
Purity	≥85% (SDS-PAGE)
Sequence	MSIDLSLPEL PILHPRITVV GVGGAGGNAV NNMIQSNLQG VNFVVANTDA QALEKSLCDK KIQLGINLTK GLGAGALPDV GKGAAEESID EIMEHIKDSH MLFITAGMGG GTGTGAAPVI AKAAREARAA VKDRAPKEKK ILTVGVVTKP FGFEGVRRMP IAELGLEELQ KYVDTLIVIP NQNLFRANE KTTFSDAFKL ADNVLHIGIR GVTDLMVMPG LINLDFADIE TVMSEMGMKAM IGTGEAEGED RAISAAEAAI SNPLLDNVSM KGAQGILINI TGGGDMTLFE VDAAANRVRE EVDENANIIF GATFDQAMEG RVRVSVLATG IDGRNNKSET SPISQSEDSE KEKFKWPYSQ SESTQDKTLE TKPAEQVSEG AKWGSNIYDI PAYLRRKK
Source	Baculovirus
Target Names	ftsZ
Protein Names	Recommended name: Cell division protein FtsZ
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.