



Recombinant 1,3-propanediol dehydrogenase (dhaT)

Product Code	CSB-YP332620DTO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P45513
Product Type	Recombinant Protein
Immunogen Species	Citrobacter freundii
Purity	>85% (SDS-PAGE)
Sequence	MSYRMFDYLV PNVNFFGPNA ISVVGERCKL LGGKKALLVT DKGLRAIKDG AVDKTLTHLR EAGIDVVVFD GVEPNPKDTN VRDGLEVFRK EHCDIIVTVG GGSPHDCGKG IGIAATHEGD LYSYAGIETL TNPLPIVAV NTTAGTASEV TRHCVLTK TKVKFVIVSW RNLPVSIND PLLMLGKPAP LTAATGMDAL THAVEAYISK DANPVTAAA IQAIRLIARN LRQAVLGSN LKARENMAYA SLLAGMAFNN ANLGYVHAMA HQLGGLYDMP HGVANAVLLP HVARYNLIAN PEKFADIAEF MGENTDGLST MDAELAIHA IARLSADIGI PQHLRDLGVK EADFPYMAEM ALKDGNAFSN PRKGNEKEIA EIFRQAF
Source	Yeast
Target Names	dhaT
Protein Names	Recommended name: 1,3-propanediol dehydrogenase EC= 1.1.1.202 Alternative name(s): 1,3-propanediol oxidoreductase 3-hydroxypropionaldehyde reductase
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.