



Recombinant Escherichia coli Uncharacterized protein yhfS (yhfS)

Product Code	CSB-BP332629ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P45545
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKTFPLQSLT IIEAQQKQFA LVDSICRHFP GSEFLTGGDL GLTPGLNQPR VTQRVEQVLA DAFHAQAAAL VQGAGTGAIK AGLAALLKPG QRLLVHDAPV YPTTRVIIIEQ MGLTLITVDF NDLSALKQVV DEQQPDAALV QHTRQQPQDS YVLADVLATL RAAGVPALTD DNYAVMKVAR IGCECGANVS TFSCFKLFGP EGVGAVVGDA DVINRIRATL YSGGSQIQGA QALEVLRGLV FAPVMHAVQA GVSERLLALL NGGAVPEVKS AVIANAQSKV LIVEFHQPIA ARVLEEAQKR GALPYPVGAE SKYEIPPLFY RLSGTFRQAN PQSEHCAIRI NPNRSGEETV LRILRESIAS I
Source	Baculovirus
Target Names	yhfS
Protein Names	Recommended name: Uncharacterized protein yhfS
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.