



Recombinant Human Ret finger protein-like 4A (RFPL4A)

Product Code	CSB-MP019602HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NLU0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAEHFKQIIR CPVCLKDLEE AVQLKCGYAC CLQCLNSLQK EPDGEGLLCR FCSVVSQKDD IKPKYKLRAL VSIIKELEPK LKSVLTMNPR MRKFQVDMTF DVDTANNYLI ISEDLRSFRS GDLSQNRKEQ AERFDTALCV LGTPRFTSGR HYWEVDVGTS QVWDVGVCKE SVNRQGKIVL SSEHGFLTVG CREGKVFAAS TVPMTPLWVS PQLHRVGIFL DVGMRSAIFY NVSDGCHIYT FIEIPVCEPW RPFFAHKRGS QDDQSILSIC SVINPSAASA PVSSEGK
Source	Mammalian cell
Target Names	RFPL4A
Protein Names	Recommended name: Ret finger protein-like 4A Alternative name(s): RING finger protein 210
Expression Region	1-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.