



Recombinant Protein parB (parB)

Product Code	CSB-MP332877ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P22997
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	MPARLGPNGY TQFDQGGGRDE GGQALPTTSI HRLQHRGAGL RAWTRRHTRP AHPRPGARGP RPARPGRVVF AGVRLAGGAG RLLQGQFRCS AAPGVPRSSA MGPRMKRRSY AMLRAAAALA VLVVASPAWA ELRGEVVRII DGDTIDLVD KQPVRVRLVD IDAPEKRQAF GERARQALAG MVFRRHVLVD EKDTDYGRGRT LGTVWVNMEL ASRPPQPRNV NAAMVHQGMA WAYRFHGAA DPEMLRLEQE ARGKRVGLWS DPHAVEPWKW RRESNNRRDE G
Source	Mammalian cell
Target Names	parB
Protein Names	Recommended name: Protein parB
Expression Region	1-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.