



# Recombinant Rice gall dwarf virus Non-structural protein 9

<b>Product Code</b>	CSB-MP332926RMJ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P23628
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rice gall dwarf virus (RGDV)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFTSSAAKTG RRRASTAKKS EDELVVDRVP NLIDGSISSG PISTGPDSQ NNKKKLPNPF APKINNLGKD GQTIEQGSML ENKSRENEGD GGSPECVTPN HTDIGTSLGQ VTTKGNECHD KYIGRGSPPN ERGGSRPSQE RGDSFVAIKY PSMLTTLAVM LNYDIHNLYE FMAVLNVTK ELNDPNVLVE RAVFDAMQYS RSRGIGDRES YSMFCYMIHG YASLMRLAEE PWSGDVSSNE SEIHIKASDM KKS VGVT LTV KPNSLWVCNK NDFARLICIF TLPDDIIAFL RTDGEDCYGG SNIYVGLDIS KPPYIPLRDV EEP
<b>Source</b>	Mammalian cell
<b>Protein Names</b>	Recommended name: Non-structural protein 9 Short name= Pns9
<b>Expression Region</b>	1-323
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.