



Recombinant Human YjeF N-terminal domain-containing protein 3 (YJEFN3)

Product Code	CSB-MP026264HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6XGL0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSSAAGPDPS EAPEERHFLR ALELQPPLAD MGRAELSSNA TTSLVQRRKQ AWGRQSWLEQ IWNAGPVCQS TAEAAALERE LLEDYRFGRQ QLVELCGHAS AVAVTKAFPL PALSRKQRTV LVVCGPEQNG AVGLVCARHL RVFEYEPTIF YPTRSLDLLH RDLTTQCEKM DIPFLSYLPT EVQLINEAYG LVVDAVLGPG VEPGEVGGPC TRALATLKLL SIPLVSLDIP SGWDAETGSD SEDGLRPDVL VSLAAPKRCA GRFSGRHHFV AGRFVPDDVR RKFALRLPGY TGTDCVAAL
Source	Mammalian cell
Target Names	YJEFN3
Protein Names	Recommended name: YjeF N-terminal domain-containing protein 3 Short name= YjeF_N3 Short name= hYjeF_N3
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.