



# Recombinant Zea mays Chalcone synthase C2 (C2)

<b>Product Code</b>	CSB-EP333027ZAX-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P24825
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Zea mays (Maize)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAGATVTVEE VRKAQRATGP ATVLAIGTAT PANCVYQADY PDYYFRITKS EHLTDLKEKF KRMCDKSMIR KRYMHLTEEF LAENPSMCAY MAPSLDARQD VVVVEVPKLG KAAAQKAIKE WGQPKSRITH LVFCTTSGVD MPGADYQLTK ALGLRPSVNR LMMYQQGCFA GGTVLRVAKD LAENNRGARV LVVCSEITAV TFRGPSESHL DSLVGQALFG DGAAAVVGA DPDDRVERPL FQLVSAAQTI LPDSEGAIDG HLREVGLTFH LLKDVPGLIS KNIGRALDDA FKPLGISDWN SIFWVAHPGG PAILDQVEAK VGLDKARMRA TRHVLSEYGN MSSACVLFIL DEMRKRSAED GQATTGEGLD WGVLFVFGPG LTVETVVLHS VPITGAATA
<b>Source</b>	E.coli
<b>Target Names</b>	C2
<b>Protein Names</b>	Recommended name: Chalcone synthase C2 EC= 2.3.1.74 Alternative name(s): Naringenin-chalcone synthase C2
<b>Expression Region</b>	1-400
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.