



Recombinant *Saccharomyces cerevisiae* Protein phosphatase 2C homolog 6 (PTC6)

Product Code	CSB-BP333093SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25646
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SRINRL ENQDNYSIKM MRSWPNAYGS ALNCSVFDGH GEKGAQLSQL LADKLCSSLD FPEPSWDKQD LKKLVQEYAR RFPEGNYWKH KLSTFEKFYN KFIKNCNSKQ ELLLMKEGDS AILGQNGGRM IFDKMGNIID KIALLTELDR LRLFYGFARF DLDQCCGLGT AAGSTASSIF LYPYDDPNAP IDEGKDDDSW IISHSGLLKL IVTQVGDSKI ILCDQDGAH ALTTTHHINS SRERHRLSID PSRLDPDAFG ETRFLNNFAN TRSFGDVAGK PYGISSEPI FSFLVGNTLH LPRSERSKLP FNGDECFLAL VTDGITNKLA DQEVVDLITS TVNSWGLKKA TPQFVAEETI KFIQAIATKH SDNATCVVVR LSNWGNWPNV DRTGPQRETK LMNAQSNETK LN
Source	Baculovirus
Target Names	PTC6
Protein Names	Recommended name: Protein phosphatase 2C homolog 6 Short name= PP2C-6 EC= 3.1.3.16
Expression Region	35-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.