



Recombinant Arabidopsis thaliana Cyclin-dependent kinase B1-1 (CDKB1-1)

Product Code	CSB-MP333114DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25859
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEKYEKLEKV GEGTYGKVYK AMEKG TGKLV ALKKTRLEMD EEGIPPTALR EISLLQMLST SIYVVRLLCV EHVHQPSTKS QSTKSNLYLV FEYLDTDLKK FIDSYRKGPN PKPLEPFLIQ KLMFQLCKGV AHCHSHGVLH RDLKPQNLLL VKDKELLKIA DLGLGRAFTV PLKSYTHEIV TLWYRAPEVL LGSTHYSTGV DMWSVGCIFA EMVRRQALFP GDSEFQQLLH IFRLLGTPTE QQWPGVSTLR DWHVYPKWEP QDLTLAVPSL SPQGVDLLTK MLKYNPAERI SAKTALDHPY FDSLDKSQF
Source	Mammalian cell
Target Names	CDKB1-1
Protein Names	Recommended name: Cyclin-dependent kinase B1-1 Short name= CDKB1;1 EC= 2.7.11.22 EC= 2.7.11.23 Alternative name(s): Cell division control protein 2 homolog B
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.