



Recombinant Salmonella typhimurium Flagellar hook-length control protein (fliK)

Product Code	CSB-EP333162SXB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P26416
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	>85% (SDS-PAGE)
Sequence	MITLPQLITT DTDMTAGLTS GKTTGSAEDF LALLAGALGA DGAQGGKDARI TLADLQAAGG KLSKELLTQH GEPGQAVKLA DLLAQKANAT DETLTDLTQA QHLLSTLTPS LKTSALAALS KTAQHDEKTP ALSDEDLASL SALFAMPLPGQ PVATPVAGET PAENHIALPS LLRGDMPSAP QEETHLSFS EHEKGKTEAS LARASDDRAT GPALTPLVVA AAATSAKVEV DSPPAPVTHG AAMPTLSSAT AQPLPVASAP VLSAPLSHE WQQTFSQQVM LFTRQGQQA QLRLHPEELG QVHISLKLDD NQAQLQMVSP HSHVRAALEA ALPMLRTQLA ESGIQLGQSS ISSEFAGQQ QSSSQQSSR AQHTDAFGAE DDIALAAPAS LQAAARGNGA VDIFA
Source	E.coli
Target Names	fliK
Protein Names	Recommended name: Flagellar hook-length control protein
Expression Region	1-405
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.