



# Recombinant Human TANK-binding kinase 1-binding protein 1 (TBKBP1)

<b>Product Code</b>	CSB-MP023233HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A7MCY6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MESMFEDDIS ILTQEALGPS EVWLDSPGDP SLGGDMCSAS HFALITAYGD IKERLGGLER ENATLRRRLK VYEIKYPLIS DFGEEHGFSL YEIKDGSLL VEKVSLQQRL NQFQHELQKN KEQEEQLGEM IQAYEKLCVE KSDLETELRE MRALVETHLR QICGLEQQLR QQQGLQDAAF SNLSPPPAPA PPCTDLDLHY LALRGGGSLG HAGWPGSTPS VSDLERRRLE EALEAAQGEA RGAQLREEQL QAECERLQGE LKQLQETRAQ DLASNQSERD MAWVKRVGDD QVNLALAYTE LTEELGRLRE LSSLQGRILR TLLQEQARSG GQRHSPLSQ R HSPAPQCPSP SPPARAAPP C PPCQSPVQR RSPVPPCPSP QRRRSPASPS C PSPVPQRRS PVPPSCQSPS PQRSPVPPS CPAPQPRPPP PPPPGERTLA ERAYAKPPSH HVKAGFQGR SYSELAEGAA YAGASPPWLQ AEAATLPKPR AYGSELYGPG RPLSPRRAFE GIRLRFKQP SEDEWAVPT SPPSPEVGTI RCASFCAGFP IPESPAATAY AHAHAQSWP SINLLMETVG SDIRSCPLCQ LGFPVGYRDD ALIKHIDSHL ENSKI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TBKBP1
<b>Protein Names</b>	Recommended name: TANK-binding kinase 1-binding protein 1 Short name= TBK1-binding protein 1
<b>Expression Region</b>	1-615
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.