



Recombinant *Pseudomonas stutzeri* P-protein (pheA)

| | |
|--------------------------|---|
| Product Code | CSB-BP333272FGJ |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P27603 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Pseudomonas stutzeri</i> (<i>Pseudomonas perfectomarina</i>) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSEADQLKAL RVRIDSLDER ILDLISERAR CAQEVARVKT ASWPKAEAEAV FYRPEREAWV LKHIMELNKG PLDNEEMARL FREIMSSCLA LEQPLRVAYL GPEGTFSQAA ALKHFGHSV I SKPMAAIDEV FREVVAGAVN FGVVPVENST EGAVNHTLDS FLEHDIVICG EVELRIHHHL LVGETTKTDR ITRIYSHAQS LAQCRKWLDA HYPNVERVAV SSNADAAKRV KSEWNSAAIA GDMAAQLYGL SKLAEKIEDR PVNSTRFLII GSQEVPTGD DKTSIIVSMR NKPGALHELL MPFHSNGIDL TRIETRPSRS GKWTYVFFID CMGHHQDPLI KNVLEKIGHE AVALKVLGSY PKAVL |
| Source | Baculovirus |
| Target Names | pheA |
| Protein Names | Recommended name: P-protein Including the following 2 domains: Chorismate mutase Short name= CM EC= 5.4.99.5 Prephenate dehydratase Short name= PDT EC= 4.2.1.51 |
| Expression Region | 1-365 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |