



Recombinant *Pseudomonas stutzeri* P-protein (pheA)

Product Code	CSB-YP333272FGJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27603
Product Type	Recombinant Protein
Immunogen Species	<i>Pseudomonas stutzeri</i> (<i>Pseudomonas perfectomarina</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSEADQLKAL RVRIDSLDER ILDLISERAR CAQEVARVKT ASWPKAEAEAV FYRPEREAWV LKHIMELNKG PLDNEEMARL FREIMSSCLA LEQPLRVAYL GPEGTFSQAA ALKHFGHSV I SKPMAAIDEV FREVVAGAVN FGVVPVENST EGAVNHTLDS FLEHDIVICG EVELRIHHHL LVGETTKTDR ITRIYSHAQS LAQCRKWLDA HYPNVERVAV SSNADAAKRV KSEWNSAAIA GDMAAQLYGL SKLAEKIEDR PVNSTRFLII GSQEVPTGD DKTSIIVSMR NKPGALHELL MPFHNSGIDL TRIETRPSRS GKWTYVFFID CMGHHQDPLI KNVLEKIGHE AVALKVLGSY PKAVL
Source	Yeast
Target Names	pheA
Protein Names	Recommended name: P-protein Including the following 2 domains: Chorismate mutase Short name= CM EC= 5.4.99.5 Prephenate dehydratase Short name= PDT EC= 4.2.1.51
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.