



# Recombinant *Oryza sativa* subsp. *japonica* Alpha-amylase isozyme 3A (AMY1.2)

<b>Product Code</b>	CSB-EP333307OFG-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P27932
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           QTQI LFQGFNWDSW KKQGGWYNML KDQVGDIASA GVTHVWLPPP            THSVSPQGYM PGRLYDLNAS KYGTKAELKS LIAAFHAKGI KCVADIVVNH            RCADDKDGRG VYCIFKGGGP RGCLDWGPSM ICCDDTQYSD            GTGHRDTGAD FAAAPDIDHL NPLVQRELSL WLRWLRRDVG FDGWRLDFAK            GYSAAVARTY VQNARPSFVV AEIWNSLSYD GDGKPAANQD            GQRQELVNWV KQVGGPATAF DFTTKGILQS AVQGELWRMR            DKDGKAPGMI GWYPEKAVTF VDNHDTGSTQ RMWPFPSDKV ILGYAYILTH            PGVPCIFYDH VFDWNLKQEI NALAATRKRN GINAGSKLRV LAAESDMYVA            MVDERVITKI GPRIDVGNII PSDFHIVAHG NDYCVWEKSG LRVPEPEGRR         </p>
<b>Source</b>	E.coli
<b>Target Names</b>	AMY1.2
<b>Protein Names</b>	Recommended name: Alpha-amylase isozyme 3A EC= 3.2.1.1 Alternative name(s): 1,4-alpha-D-glucan glucanohydrolase
<b>Expression Region</b>	27-440
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.